



Mentype® **DIPscreen** **Manual**

The entry to a quantitative chimerism analysis

In-Vitro-Diagnostics



25
100
400



Version DISGAv2en



45-45410-0025
45-45410-0100
45-45410-0400



Batch Code



Biotype Diagnostic GmbH
Moritzburger Weg 67
D-01109 Dresden
Germany

Made in Germany

Biotype Diagnostic GmbH develops, produces and markets their PCR-based rapid Mentype® Detection Kits. Our products provide customers with fast and reliable testing methods for professional medical diagnostics.

Our Mentype® Test Kits guarantee highest quality standards for clinical research and diagnostics.

For information and enquiries about the Mentype® **DIPscreen** please do not hesitate to get in touch or visit www.biotype.de

Product description

Mentype® **DIPscreen** is a multiplex-PCR application developed to identify DIP polymorphisms that individually occur in donor or recipient, and, constitute informative loci. In a single multiplex-PCR 33 DIP-loci are simultaneously screened together with the gender specific locus Amelogenin. Mentype® **DIPscreen** is a multiplex-PCR application that mediates monitoring of chimerism samples after stem cell transplantation. The flexible assay format allows an individual diagnostics at any time required.

Analysis of molecular chimerism resulting from allogeneic stem cell transplantation has become a well established method to control the course of transplant engraftment and to assess the risk of threatening relapse. Molecular chimerism analysis can be performed on diverse DNA-sequence motifs of which biallelic short insertion/deletion polymorphisms (DIPs, INDELs) offer substantial benefits. Polymerase-mediated amplification of DIP-markers does not result in formation of stutter peaks that can hamper clear analysis. Moreover, these polymorphisms are best suited for allele specific quantitative approaches. Mentype® **DIPscreen** is a DIP-based chimerism analysis and therefore accounts for an unambiguous donor/recipient differentiation and highly clear chimerisms monitoring.

The 33 DIP loci addressed by Mentype® **DIPscreen** are distributed over 18 chromosomes, and are at least separated by 10 Mbp each (see Tab. 1). The detection limit of Mentype® **DIPscreen** is about **200 pg genomic DNA**. The optimal range under standard conditions is **1.0 -2.0 ng DNA**. For fast and sensitive fragment length analysis primers are fluorescence-labelled with **6-FAM, BTG, or BTY**.

The test kit was validated and evaluated using the GeneAmp® 9700 Silver, Eppendorf Mastercycler ep-S, Biometra T1, ABI PRISM® 3130 Genetic Analyzer running with 36 cm capillary array and POP4® polymer. Development, manufacture and distribution of Biotype® products are certified according to DIN EN ISO13485.

Content

| | |
|--|-----------|
| 1. Description of the Mentype® DIPscreen | 6 |
| 2. PCR amplification | 10 |
| 2.1 Master mix preparation | 10 |
| 2.2 PCR amplification parameter | 11 |
| 3. Electrophoresis using the ABI PRISM® 310 Genetic Analyzer | 12 |
| 3.1 Matrix generation | 12 |
| 3.2 Sample preparation | 15 |
| 3.3 Setting up the Data Collection Software | 15 |
| 3.4 Analysis parameter | 16 |
| 4. Electrophoresis using the ABI PRISM® 3100-Avant/ 3100 Genetic Analyzer | 17 |
| 4.1 Spectral calibration / matrix generation | 17 |
| 4.2 Sample preparation | 19 |
| 4.3 Setting up the Data Collection Software | 20 |
| 4.4 Analysis parameter / analysis method | 21 |
| 5. Electrophoresis using the ABI PRISM® 3130/ 3130xl Genetic Analyzer | 22 |
| 5.1 Spectral calibration / matrix generation | 22 |
| Sample preparation | 25 |
| 5.2 Setting up the Data Collection Software | 26 |
| 5.3 Analysis parameter / analysis method | 28 |
| 6. Electrophoresis using the ABI PRISM® 3500/ 3500xL Genetic Analyzer | 29 |
| 6.1 Spectral calibration / matrix generation | 29 |
| 6.2 Sample preparation | 32 |
| 6.3 Setting up a run | 33 |
| 7. Analysis | 36 |
| 7.1 Biotype® template files | 37 |
| 7.2 Controls | 38 |
| 7.3 Lengths of fragments and alleles | 39 |
| 8. Interpretation of results | 43 |
| 9. References | 44 |
| 10. Explanation of Symbols | 45 |

| | |
|---|-----------|
| A Analytical Validation | 46 |
| A a) Determination of the Standard Reaction and batch-specific Tolerance..... | 46 |
| A b) Genotyping Accuracy | 46 |
| A c) Analytical Specificity..... | 47 |
| A d) Analytical Sensitivity..... | 47 |
| A e) Assays Performance with Different PCR-Thermocyclers..... | 47 |
| A f) Mixed DNA Samples | 48 |
| A g) PCR Annealing Temperatures | 48 |
| A h) Fluctuation of PCR buffer Batches | 48 |
| A i) In-use Stability..... | 49 |
| B Clinical Performance Data | 49 |
| B a) Study Design, ethics and regulatory aspects | 49 |
| B b) Reference Methods..... | 49 |
| B c) DNA-Extraction and Purification | 49 |
| B d) Results..... | 50 |
| B e) References | 52 |

1. Description of the Mentype® DIPscreen

Table 1. Locus-specific information of Mentype® DIPscreen

| DIP Locus | Chromosomal position | Motive (-DIP / +DIP) |
|------------------|----------------------|--------------------------------------|
| FAM Panel | | |
| AM X | Xp22.1-22.3 | |
| AM Y | Yp11.2 | |
| HLD106 | 16q13 | -/AATGCGT |
| HLD70 | 6q16.1 | -/AGCA |
| HLD84 | 8q24.12 | -/CTTTC |
| HLD103 | 12q23.1 | -/GCTTATAA |
| HLD104 | 13q32.1 | -/ACTC |
| HLD116 | 18p11.22 | -/AGGTGTGGAACAACATGATAC |
| HLD112 | 17p12 | -/TTGTA |
| HLD307 | Xp11.23 | -/TCAACCAA |
| HLD310 | 2p22.3 | -/GTCTGGTT |
| HLD110 | 16q22.1 | -/TCCCTG |
| HLD133 | 3p22.1 | -/CAACCTGGATT |
| HLD79 | 7q31.2 | -/AATCT |
| HLD105 | 14q24.3 | -/ATAGACAA |
| HLD140 | 3q23 | -/GGTAGTATGGGCCCT |
| HLD163 | 12q24.31 | -/AACTACGGCACGCC |
| BTG Panel | | |
| HLD91 | 11q14.1 | -/GATA |
| HLD23 | 18p11.32 | -/CTTTAA |
| HLD88 | 9q22.33 | -/CCACAAAGA |
| HLD101 | 15q26.1 | -/GTAG |
| HLD67 | 5q33.3 | -/CTACTGAC |
| HLD301 | 17q21.32 | -/CAGGGGCTC |
| HLD53 | 3q22.1 | -/ATGT |
| HLD97 | 13q13.1 | -/AGAGAAAGCTGAAG |
| HLD152 | 16p13.2 | -/TGGTCAAAGGCA |
| HLD128 | 1q31.3 | -/ATTAATA |
| HLD134 | 5q11.2 | -/ATGATGGTTCTTCAGA |
| HLD305 | 20q11.22 | -/CAAGGTCCCACCACACTCGCGTGGGA |
| BTY Panel | | |
| HLD48 | 2q11.2 | -/GACTT |
| HLD114 | 17p13.2 | -/TCCTATTCTACTCTGAAT |
| HLD304 | 9q34.3 | -/GAGCTGCTCAAGAGAGAGG |
| HLD131 | 7q36.2 | -/TTGGGCTTATT |
| HLD38 | 1q32.2 | -/TAGTT |
| HLD82 | 7q21.3 | -/ACCTCTACTCCTGGTCTATTCTGGTCACATGACT |

Abbreviations: HLD = Human Locus DIP, -DIP = Deletion, +DIP = Insertion

Table 1 shows the chromosomal position, motif and respective reference allele of DIP-loci addressed by Mentype® DIPscreen.

Kit Content

Mentype® DIPscreen (100 Reactions)

| | |
|-----------------------------|--------|
| Nuclease-free water | 3.0 ml |
| Reaction mix A | 500 µl |
| Primer mix | 500 µl |
| Multi Taq2 DNA polymerase | 60 µl |
| Control DNA XY13 (2ng/µl) | 10 µl |
| DNA Size Standard 550 (BT0) | 50 µl |
| Allelic ladder | 25 µl |

Ordering information

| Product | Reactions | Order number |
|--------------------|---------------|---------------|
| Mentype® DIPscreen | 25 reactions | 45-45410-0025 |
| Mentype® DIPscreen | 100 reactions | 45-45410-0100 |
| Mentype® DIPscreen | 400 reactions | 45-45410-0400 |

Storage

Store all components at -20 °C and avoid repeated thawing and freezing. Primer mix and allelic ladder must be stored protected from light. The DNA samples and post-PCR reagents (allelic ladder and DNA size standard) should be stored separately from the PCR reagents. The expiry date is indicated on the kit cover.

Additionally required reagents

Additional reagents are required in order to use the Biotype® PCR Amplification Kit:

| Reagent | Supplier | Order Number |
|--|-------------------------|---------------|
| Hi-Di™ Formamide, 25 ml | Applied Biosystems | 4311320 |
| Matrix Standards BT5 single-capillary instruments (5x25 µl) | Biotype Diagnostic GmbH | 00-10411-0025 |
| Matrix Standards BT5 multi-capillary instruments (25 µl) | Biotype Diagnostic GmbH | 00-10421-0025 |
| Matrix Standards BT5 multi-capillary instruments (50 µl) | Biotype Diagnostic GmbH | 00-10421-0050 |

Warning and safety instructions

The PCR Amplification Kit contains the following potentially hazardous chemicals:

| Kit component | Chemical | Hazards |
|----------------------|-----------------------------|---|
| Reaction mix | Sodium azide NaN_3 | toxic if swallowed, develops toxic gases when it gets in contact with acids |


Observe the Material Safety Data Sheets (MSDS) for all Biotype® products, which are available on request. Please contact the respective manufacturers for copies of the MSDS for any additionally needed reagents.

Quality assurance

All kit components undergo an intensive quality assurance process at Biotype Diagnostic GmbH. The quality of the test kits is permanently monitored to ensure unrestricted usability. Please contact us if you have any questions regarding quality assurance.

Trademarks and patents

Mentype® is a registered trademark of Biotype Diagnostic GmbH.
 ABI PRISM®, GeneMapper®, GeneAmp® and Applied Biosystems® are registered trademarks of Applied Biosystems LLC.
 Under the law of Europe POP4® is registered trademark of Applied Biosystems LLC.
 The PCR is covered by patents. Patentees are Hoffmann-La Roche Inc. and F. Hoffmann-La Roche (Roche).

Outline of working steps performed with Mentype® DIP-products

From sample to analysis – Monitoring chimerism with the Mentype® DIPscreen

Protocols for PCR amplification, electrophoresis and analysis

2. PCR amplification

2.1 Master mix preparation

The table below shows the volumes of all PCR reagents per 25 µl reaction volume, including a sample volume of 1.0 µl (template DNA). The number of reactions to be set up shall be determined taking into account positive and negative control reactions. Add one or two reactions to this number to compensate the pipetting error.

| Component | Volume |
|---|---------|
| Nuclease-free water | 13.4 µl |
| Reaction mix A* | 5.0 µl |
| Primer mix | 5.0 µl |
| Multi Taq2 DNA polymerase (hot start, 2.5 U/µl) | 0.6 µl |
| Volume of master mix | 24.0 µl |

* contains Mg²⁺, dNTPs, BSA

All components should be mixed (vortex) and centrifuged for about 10 s before preparing the master mix. The DNA volume applied to the assay depends on its concentration. For reference samples 1 µl is mostly sufficient. For critical patient samples the amount of template can be increased appropriately. Fill up the final reaction volume to 25 µl with nuclease-free water.

Generally, DNA templates shall be stored in nuclease-free water or in diluted TE buffer (10 mM Tris HCl, pH 8.0 and 1 mM EDTA), e.g. 0.1x TE buffer.

The primer mixes are adjusted for balanced peak heights at **28 PCR cycles** and **1 ng Control DNA XY13** in a reaction volume of 25 µl. If more DNA template is applied, higher peaks can be expected for small PCR fragments and relatively low peaks for large fragments. Reduce the amount of DNA template to correct this imbalance.

Positive control

For the positive amplification control, dilute the Control DNA XY13 to 1 ng/µl. Instead of the template DNA, pipette the diluted Control DNA into a reaction tube containing the PCR master mix.

Negative control

For the negative amplification control, pipette nuclease-free water instead of template DNA into a reaction tube that contains the PCR master mix.

Template DNA

Sometimes, measured DNA concentration varies depending on the quantification method used. It might thus be necessary to adjust the optimal DNA amount.

2.2 PCR amplification parameter

Perform a “hot start” PCR in order to activate the Multi Taq2 DNA Polymerase and to prevent formation of non-specific amplification products.

The number of cycles depends on the amount of DNA applied. 28 PCR cycles are recommended for all samples.

Standard method

Recommended for all DNA samples

| Temperature | Time | |
|-------------|---|------------------|
| 94°C | 4 min (hot start for activation of the Multi Taq2 DNA polymerase) | |
| 94°C | 30 s | |
| 60°C | 120 s | 28 cycles |
| 72°C | 75 s | |
| 68°C | 60 min | |
| 10°C | ∞ | hold |

Note: To provide an optimal kit balance the ramping rate of the thermal cycler should be adjusted to 4-5 °C/s.

Very small amounts of DNA may result in statistical dropouts and imbalances of the peaks. Increasing numbers of PCR cycles raise the risk of cross contamination caused by minimal amounts of impurities. Furthermore, unspecific amplification products could appear.

3. Electrophoresis using the ABI PRISM® 310 Genetic Analyzer

For general instructions on instrument setup, matrix generation and application of the GeneScan® or GeneMapper® ID software, refer to the *ABI PRISM® 310 Genetic Analyzer User's Manual*. Electrophoresis using the GeneScan® software is described below.

The virtual **filter set G5** shall be used for combined application of the five fluorescent labels **6-FAM, BTG, BTY, BTR, and BTO** (the matrix standard will be called **BT5** hereinafter).

Material

| | |
|-----------|---------------------------------------|
| Capillary | 47 cm / 50 µm (green) |
| Polymer | POP4® for 310 Genetic Analyzer |
| Buffer | 10x Genetic Analyzer Buffer with EDTA |

3.1 Matrix generation

Prior to conducting DNA fragment size analysis with the **Filter Set G5**, a matrix with the five fluorescent labels **6-FAM, BTG, BTY, BTR, and BTO** must be generated.

| Color | Matrix standard |
|------------|-----------------|
| Blue (B) | 6-FAM |
| Green (G) | BTG |
| Yellow (Y) | BTY |
| Red (R) | BTR |
| Orange (O) | BTO |

Five electrophoresis runs shall be conducted, one for each fluorescent label, **6-FAM, BTG, BTY, BTR, and BTO**. Use the same conditions as for samples and allelic ladders of the Biotype® test kit to generate suitable matrix files.

| Matrix sample | Component | Volume |
|--|------------------------------|---------|
| Matrix sample 1 | Hi-Di™ Formamide | 12.0 µl |
| | Matrix standard 6-FAM | 1.0 µl |
| Matrix sample 2 | Hi-Di™ Formamide | 12.0 µl |
| | Matrix standard BTG | 1.0 µl |
| Matrix sample 3 | Hi-Di™ Formamide | 12.0 µl |
| | Matrix standard BTY | 1.0 µl |
| Matrix sample 4 | Hi-Di™ Formamide | 12.0 µl |
| | Matrix standard BTR | 1.0 µl |
| Matrix sample 5 | Hi-Di™ Formamide | 12.0 µl |
| | Matrix standard BTO | 1.0 µl |
| - Denaturation for 3 min at 95°C | | |
| - Cool down to 4°C and place samples on the autosampler tray | | |

- Create a **Sample Sheet** choose **5 Dyes** and enter a sample designation

Injection list for matrix generation

| Parameter | Set up |
|----------------|------------------------------|
| Module File | GS STR POP4 (1 ml) G5 |
| Matrix File | NONE |
| Size Standard* | NONE |
| Injection [s] | 5 |
| Injection [kV] | 15.0 |
| Run [kV] | 15.0 |
| Run [°C] | 60 |
| Run Time [min] | 24 |

* Prepare matrix standards always **without DNA Size Standard (BTO)**

Analysis of the matrix samples

- Run the GeneScan[®] software
- **File** → **New** → **Project** (open folder of current run) → **Add Sample Files**
- Select a matrix sample in the **Sample File** column
- **Sample** → **Raw Data**
- Check the matrix samples for a flat baseline. As shown in the figure below there should be at least five peaks with peak heights about 1000-4000 RFU (Y-axis) for each matrix sample (optimal range: 2000-4000 RFU)




Fig. 1 Electropherogram with raw data of the matrix standard 6-FAM

- Select an analysis range with flat baseline and re-inject the matrix sample if necessary
- Note down start and end value (data points) of the analysis range, e.g. start value 3200, end value 5500
- Calculate the difference, e.g. $5500 - 3200 = 2300$ data points

Generation of a new matrix

- **File → New → Matrix**




Fig. 2 Matrix sample selection

- Import matrix samples for all dyes (B, G, Y, R, O)
- Enter a **Start At** value, e.g. 3200
- Enter the calculated difference under **Points**, e.g. 2300
- Click on **OK** to calculate the new matrix

| | B | G | Y | R | O |
|---|--------|--------|--------|--------|--------|
| B | 1.0000 | 0.1811 | 0.0051 | 0.0418 | 0.0006 |
| G | 0.6891 | 1.0000 | 0.2056 | 0.3259 | 0.0017 |
| Y | 0.4687 | 0.8068 | 1.0000 | 0.9119 | 0.0029 |
| H | 0.1944 | 0.3619 | 0.5311 | 1.0000 | 0.0095 |
| O | 0.0160 | 0.0304 | 0.0477 | 0.2082 | 1.0000 |

Fig. 3 New matrix BT5

- Save the matrix in the matrix folder: **File → Save as**, e.g. Matrix BT5

Matrix check

Check the new matrix with current samples.

- **File → New → Project** (open folder of the respective run) → **Add Sample Files**
- Select sample(s) in the **Sample File** column
- **Sample → Install New Matrix** (open matrix folder and select new matrix)
- Re-analyse your samples

There should be no pull-up peaks between the dye panels (B, G, Y, R, O) with the new matrix.

3.2 Sample preparation

| Component | Volume |
|--|---------|
| Hi-Di™ Formamide | 12.0 µl |
| DNA Size Standard 550 (BTO) | 0.5 µl |
| Prepare 12 µl of the mix (formamide + DNA size standard) for all samples | |
| Add 1 µl PCR product (diluted if necessary) or allelic ladder | |
| - Denaturation for 3 min at 95 °C | |
| - Cool down to 4 °C and place samples on the autosampler tray | |

Signal intensities

Options to increase the signal intensity:

- Reduce the volume of the DNA Size Standard 550 (BTO) to peak heights of about 500 relative fluorescent units (RFU)
- Purify the PCR products before starting the analysis

3.3 Setting up the Data Collection Software

- Create a **Sample Sheet** and enter a sample designation

Injection list

| Parameter | Set up |
|------------------|------------------------------|
| Module File | GS STR POP4 (1 ml) G5 |
| Matrix File | e.g. Matrix BT5 |
| Size Standard | e.g. SST-BTO_60-450bp |
| Injection [s]* | 5 |
| Injection [kV] | 15.0 |
| Run [kV] | 15.0 |
| Run [°C] | 60 |
| Run Time [min]** | 26 |

* Deviating from the standard settings, the injection time may range between 1 and 20 s depending on the type of sample. If reference samples with very high signal intensities are recorded, a shorter injection time may be selected in order to avoid pull-up peaks. For samples with low DNA content an injection time of up to 20 s may be necessary.

** Depending on the analysis conditions, the run time for Mentype® **DIPscreen** was modified in order to be able to analyse fragments with lengths of up to **450 bp**.

3.4 Analysis parameter

The recommended analysis parameters are:

| | |
|-----------------------|---|
| Analysis Range | Full Range |
| Data Processing | Baseline: Checked Multicomponent: Checked Smooth Options: Light |
| Peak Detection | Peak Amplitude Thresholds B:* Y:* G:* R:* O:* Min. Peak Half Width: 2 pts Polynomial Degree: 3 Peak Window Size: 15 pts |
| Size Call Range | Min: 60 Max: 550 |
| Size Calling Method | Local Southern Method |
| Split Peak Correction | None |

* The peak amplitude threshold (cutoff value) corresponds to the minimum peak height that will be detected by the GeneScan[®] or GeneMapper[®] ID software. Thresholds are usually 50-200 RFU and should be determined individually by the laboratory. Recommendation: The minimal peak height should be three times as high as the background noise of the baseline.

4. Electrophoresis using the ABI PRISM® 3100-Avant/3100 Genetic Analyzer

For detailed instructions on instrument setup, spectral calibration, application of the ABI PRISM® 3100 Data Collection Software version 1.01 or 1.1 and the GeneScan® software, refer to the *ABI PRISM® 3100-Avant/3100 Genetic Analyzer User's Manual*. For systems with Data Collection Software 2.0 or 3.0 refer to chapter 5.

The system with 4 capillaries is named ABI 3100-Avant, and the system with 16 capillaries is named ABI 3100.

The virtual **filter set G5** shall be used for combined application of the five fluorescent labels **6-FAM, BTG, BTY, BTR, and BTO** (the matrix standard will be called **BT5** hereinafter).

Material

| | |
|-----------|---|
| Capillary | 36 cm Capillary Array for 3100-Avant/3100 |
| Polymer | POP-4® Polymer for 3100 |
| Buffer | 10x Genetic Analyzer Buffer with EDTA |

4.1 Spectral calibration / matrix generation

Proper spectral calibration is critical to evaluate multicolour systems with the ABI PRISM® 3100-Avant/3100 Genetic Analyzer and shall be done prior to conducting fragment length analysis with the five fluorescent labels **6-FAM, BTG, BTY, BTR, and BTO**. The calibration procedure creates a matrix which is used to correct the overlapping of fluorescence emission spectra of the dyes.

Spectral calibration comprises the following steps:

- Preparation of the spectral calibration standards
- Loading the standards to the 96-well reaction plate (one sample per capillary)
- Entering the plate composition
- Performing a spectral calibration run and checking the matrix

Setting up the spectral calibration standards

Example for 4 capillaries / ABI 3100-Avant

| Component | Volume |
|--|---------|
| Hi-Di™ Formamide | 60.0 µl |
| Matrix standard BT5 | 5.0 µl |
| - Load 12 µl of the mix to a 96-well reaction plate, e.g. position A1–D1 | |
| - Denaturation for 3 min at 95 °C | |
| - Cool down to 4 °C and place samples on the autosampler tray | |

Example for 16 capillaries / ABI 3100

| Component | Volume |
|--|----------|
| Hi-Di™ Formamide | 204.0 µl |
| Matrix standard BT5 | 17.0 µl |
| - Load 12 µl of the mix to a 96-well reaction plate, e.g. position A1-H1 and A2-H2 | |
| - Denaturation for 3 min at 95 °C | |
| - Cool down to 4 °C and place samples on the autosampler tray | |

Performing a spectral calibration run

First of all, the parameter file for **DyeSetG5** must be modified once to achieve successful calibration with the Data Collection software version 1.0.1 or 1.1.

Spectral parameter

To change settings in the parameter file go to the following path:

D:\AppliedBio\Support Files\Data Collection Support Files\CalibrationData\Spectral Calibration\ParamFiles

- Select **MtxStd(Genescan_SetG5)** to open the PAR-file
- Change **Condition Bounds Range** to [1.0; 20.0]
- Select **File → Save As** to save the parameter file under a new name, e.g. MtxStd(Genescan_SetG5_BT5).par

Always use this parameter file for spectral calibration runs using Biotype® matrix standards **BT5**.

Plate Editor for spectral calibration (I)

- Place the 96-well plate on the autosampler tray
- Run the ABI PRISM® 3100 Data Collection software
- In **Plate View** click **New** to open the **Plate Editor** dialog box
- Enter a name of the plate
- Select **Spectral Calibration**
- Select **96-Well** as plate type and click on **Finish**

Plate editor for spectral calibration (II)

| Parameter | Set up |
|---------------------|--|
| Sample Name | Type name for the matrix samples |
| Dye Set | G5 |
| Spectral Run Module | <i>Default</i> (e.g. Spect36_POP4®) |
| Spectral Parameters | MtxStd(GeneScan_SetG5_BT5).par (parameters created before) |

- Click into the column header to select the entire column, select **Edit** → **Fill Down** to apply the information of the selected samples and confirm with **OK**
- Link your reaction plate on the autosampler tray with the created plate ID and start run
- On completion of the run check in the **Spectral Calibration Result** dialog box if all capillaries have successfully passed calibration (label **A**). If individual capillaries are labelled **X**, refer to *ABI PRISM® Genetic Analyzer User's Manual*.
- Click on **OK** to confirm completion of the run

Matrix check

- Select **Tools** → **Display Spectral Calibration** → **Dye Set** → **G5** to review the spectral calibration profile for each capillary
- The quality value (**Q value**) must be greater than 0.95 and the condition number (**C value**) must be between 1 and 20. Both values must be within the previously determined range
- Check the matrix samples for a flat baseline. There should be five peaks with peak heights of about 1000-5000 RFU (Y-axis) in each matrix sample (optimal range: 2000-4000 RFU)
- If all capillaries have passed the calibration, the last calibration file for **Dye Set G5** must be activated manually under **Tools** → **Set Active Spectral Calibration**. Rename the calibration file under **Set Matrix Name** (e.g. BT5_Date of calibration)
- If calibration was not successful, try to re-inject the samples with higher injection voltage or injection time. The editing of the Spectral Run Module will be necessary. You can re-inject the same samples up to three times. Otherwise use more matrix standard for spectral calibration
- Check the new matrix with your current samples. There should be no pull-up peaks between the dye panels (B, G, Y, R, O) with the new matrix

4.2 Sample preparation

| Component | Volume |
|--|---------|
| Hi-Di™ Formamide | 12.0 µl |
| DNA Size Standard 550 (BTO) | 0.5 µl |
| Prepare 12 µl of the mix (formamide + DNA size standard) for all samples | |
| Add 1 µl PCR product (diluted if necessary) or allelic ladder | |
| - Denaturation for 3 min at 95 °C | |
| - Cool down to 4 °C and place the samples on the autosampler tray | |

Since injections take place simultaneously on all capillaries, 4 or 16 samples must be pipetted on the plate of multi-capillary analyzers. If fewer samples are analysed, the empty positions must be filled with 12 µl Hi-Di™ Formamide.

To ensure a reliable allelic assignment on multi-capillary analyzers, several ladders should be run.

Room temperature may influence the performance of PCR products on multi-capillary instruments, so that shoulder peaks or split peaks occur especially at low temperatures. Pay attention to keeping ambient conditions as recommended by the instrument manufacturer. Optimal will be a stable room temperature > 22 °C.

Signal intensities

Options to increase the signal intensity:

- Reduce the volume of the DNA Size Standard 550 (BTO) to peak heights of about 500 relative fluorescent units (RFU)
- Purify the PCR products before starting the analysis

4.3 Setting up the Data Collection Software

Edit the default run module in **Dye Set G5** once for the first run.

- Select **Module Editor** to open the dialog box
- Select the appropriate **Run Module** as template from the **GeneScan** table
- Modify the **Injection Voltage** to 3 kV and the **Injection Time** to 10 s

Run Module 3kV_10s_450bp

| Parameter | Set up |
|------------------------|----------------|
| Run Temperature [°C] | <i>Default</i> |
| Cap Fill Volume | <i>Default</i> |
| Maximum Current [A] | <i>Default</i> |
| Current Tolerance [A] | <i>Default</i> |
| Run Current [A] | <i>Default</i> |
| Voltage Tolerance [kV] | <i>Default</i> |
| Pre Run Voltage [kV] | <i>Default</i> |
| Pre Run Time [s] | <i>Default</i> |
| Injection Voltage [kV] | 3.0 |
| Injection Time [s]* | 10 |
| Run Voltage [kV] | <i>Default</i> |
| Number of Steps | <i>Default</i> |
| Voltage Step Interval | <i>Default</i> |
| Data Delay Time [s] | <i>Default</i> |
| Run Time [min]** | 25 |

* Deviating from the standard settings, the injection time may range between 1 and 20 s depending on the type of sample. If reference samples with very high signal intensities are recorded, a shorter injection time may be selected in order to avoid pull-up peaks. For samples with low DNA content an injection time of up to 20 s may be necessary.

** Depending on the analysis conditions, the Run Time for Mentype® **DIPscreen** was modified in order to be able to analyse fragments with lengths of up to **450 bp**.

- Click on **Save As**, enter the name of the new module (e.g. 3kV_10s_450bp) and confirm with **OK**
- Click on **Close** to exit the **Run Module Editor**

Starting the run

- Place the prepared 96-well plate on the autosampler tray
- Run the ABI PRISM® 3100 Data Collection software
- In **Plate View** click on **New** to open the **Plate Editor** dialog box
- Enter a name of the plate
- Select **GeneScan**
- Select **96-Well** as plate type and click on **Finish**

Plate Editor

| Parameter | Set up |
|-------------------|---------------------|
| Sample Name | enter a name |
| Dyes | 0 |
| Colour Info | Ladder or sample |
| Project Name | e.g. 3100_Project1 |
| Dye Set | G5 |
| Run Module* | 3kV_10s_450bp |
| Analysis Module 1 | DefaultAnalysis.gsp |

* parameter see above

- Complete the table in the **Plate Editor** and click on **OK**
- Click into the column header to select the entire column and select **Edit → Fill Down** to apply the information of the selected samples
- Link your reaction plate on the autosampler tray with the created plate ID and start the run
- On completion of the run, view data as **Color Data** in **Array View** of the 3100 Data Collection software or as **Analyzed Sample Files** under
D:/AppliedBio/3100/DataExtractor/ExtractRuns

4.4 Analysis parameter / analysis method

The recommended analysis parameters are:

| | |
|-----------------------|---|
| Analysis Range | Full Range |
| Data Processing | Baseline: Checked Multicomponent: Checked Smooth Options: Light |
| Peak Detection | Peak Amplitude Thresholds B:* Y:* G:* R:* O:* Min. Peak Half Width: 2 pts Polynomial Degree: 3 Peak Window Size: 15 pts |
| Size Call Range | Min: 60 Max: 550 |
| Size Calling Method | Local Southern Method |
| Split Peak Correction | None |

* The peak amplitude threshold (cutoff value) corresponds to the minimum peak height that will be detected by the GeneScan® or GeneMapper® ID software. Thresholds are usually 50-200 RFU and should be determined individually by the laboratory. Recommendation: The minimal peak height should be three times as high as the background noise of the baseline.

5. Electrophoresis using the ABI PRISM® 3130/3130xl Genetic Analyzer

For detailed instructions on instrument setup, spectral calibration, or application of the ABI PRISM® Data Collection software version 3.0 and the GeneMapper® ID/ID-X software, refer to the *ABI PRISM® 3130/3130xl Genetic Analyzers Getting Started Guide*.

The system with 4 capillaries is named ABI 3130, and the system with 16 capillaries is named ABI 3130xl.

The virtual **filter set Any5Dye** shall be used for the combined application of the five fluorescent labels **6-FAM, BTG, BTY, BTR, and BTO** (the matrix standard will be called **BT5** hereinafter).

Material

| | |
|-----------|---------------------------------------|
| Capillary | 36 cm Capillary Array for 3130/3130xl |
| Polymer | POP4® Polymer for 3130 |
| Buffer | 10x Genetic Analyzer Buffer with EDTA |

5.1 Spectral calibration / matrix generation

Prior to conducting DNA fragment size analysis, it is necessary to perform a spectral calibration with the five fluorescent labels **6-FAM, BTG, BTY, BTR, and BTO** for each analyzer. The calibration procedure creates a matrix which is used to correct the overlapping of fluorescence emission spectra of the dyes.

Spectral calibration comprises the following steps:

- Preparation the spectral calibration standards
- Loading the standards to the 96-well reaction plate (one sample per capillary)
- Creating the instrument protocol for spectral calibration (Protocol Manager)
- Defining the plate composition in the plate editor (Plate Manager)
- Performing a spectral calibration run and checking the matrix

Setting up the spectral calibration standards

Example for 4 capillaries/ABI 3130

| Component | Volume |
|---|---------|
| Hi-Di™ Formamide | 60.0 µl |
| Matrix standard BT5 | 5.0 µl |
| <ul style="list-style-type: none"> - Load 12 µl of the mix to a 96-well reaction plate, e.g. position A1-D1 - Denaturation for 3 min at 95 °C - Cool down to 4 °C and place samples in the autosampler tray | |

Example for 16 capillaries/ABI 3130xl

| Component | Volume |
|--|----------|
| Hi-Di™ Formamide | 204.0 µl |
| Matrix standard BT5 | 17.0 µl |
| <ul style="list-style-type: none"> - Load 12 µl of the mix to a 96-well reaction plate, e.g. position A1-H1 and A2-H2 - Denaturation for 3 min at 95 °C - Cool down to 4 °C and place samples in the autosampler tray | |

Performing a spectral calibration run

- Place the 96-well plate on the autosampler tray
- In the **Protocol Manager** of the Data Collection software click on **New in Instrument Protocol** to open the **Protocol Editor** dialog box

Instrument Protocol for spectral calibration

| Protocol Editor | Set up |
|-----------------|--|
| Name | <i>User</i> (e.g. Spectral36_POP4_BT5) |
| Type | SPECTRAL |
| Dye Set | Any5Dye |
| Polymer* | <i>User</i> (e.g. POP4) |
| Array Length* | <i>User</i> (e.g. 36cm) |
| Chemistry | Matrix Standard |
| Run Module* | <i>Default</i> (e.g. Spect36_POP4_1) |

* Depends on the type of polymer and length of capillary used

- Click on **OK** to leave the **Protocol Editor** dialog box
- In the **Plate Manager** of the Data Collection software, click on **New** to open the **New Plate Dialog** box

Plate Editor for spectral calibration (I)

| New Plate Dialog | Set up |
|----------------------------|------------------------|
| Name | e.g. Spectral_BT5_date |
| Application | Spectral Calibration |
| Plate Type | 96-Well |
| Owner Name / Operator Name | ... |

- Click on **OK**. A new table in the **Plate Editor** will open automatically

Plate Editor for spectral calibration (II)

| Parameter | Set up |
|-----------------------|--|
| Sample Name | Enter name for the matrix samples |
| Priority | e.g. 100 |
| Instrument Protocol 1 | Spectral36_POP4_BT5 (setting described before) |

- Click into the column header to select the entire column, select **Edit** → **Fill Down** to apply the information to all selected samples, and click on **OK**
- In the **Run Scheduler** click on **Find All**, select **Link** to link the reaction plate on the autosampler up with the newly created plate record (position A or B) and start the run

GA Instruments - ga3130 - 3130-1 - Spectral Viewer




Fig. 4 Electropherogram of spectral calibration with matrix standard BT5 on an ABI 3130

Matrix check

- The quality value (**Q value**) of each capillary must be greater than 0.95 and the condition number range (**C value**) must be between 1 and 20
- Check the matrix samples for a flat baseline. As shown in the figure above, there should be five peaks with peak heights of about 1000-5000 RFU (Y-axis) in each matrix sample (optimal range: 2000-4000 RFU)
- If all capillaries have passed the test, the last calibration file for the Dye Set **Any5Dye** is activated automatically in the **Spectral Viewer**. **Rename** the calibration file (e.g. BT5_Date of calibration) using the respective button
- If calibration was not successful, try to re-inject the samples with higher injection voltage or injection time. Editing of the Spectral Run Module will be necessary. You can re-inject the same samples up to three times. Otherwise use more matrix standard for spectral calibration
- Check the new matrix with your current samples. There should be no pull-up peaks between the dye panels (B, G, Y, R, O) with the new matrix

Sample preparation

| Component | Volume |
|--|---------|
| Hi-Di™ Formamide | 12.0 µl |
| DNA Size Standard 550 (BTO) | 0.5 µl |
| Prepare 12 µl of the mix (formamide + DNA size standard) for all samples | |
| Add 1 µl PCR product (diluted if necessary) or allelic ladder | |
| - Denaturation for 3 min at 95 °C | |
| - Cool down to 4 °C and place the samples on the tray | |

Since injections take place simultaneously on all capillaries, 4 or 16 samples must be pipetted on the plate of multi-capillary analyzers. If fewer samples are analysed, the empty positions must be filled with 12 µl Hi-Di™ Formamide.

To ensure a reliable allelic assignment on multi-capillary analysers, several ladders should be run.

Room temperature may influence the performance of PCR products on multi-capillary instruments, so that shoulder peaks or split peaks occur especially at low temperatures. Pay attention to keeping ambient conditions as recommended by the instrument manufacturer. Optimal will be a stable room temperature > 22 °C.

Signal intensities

Options to increase the signal intensity:

- Reduce the volume of the DNA Size Standard 550 (BTO) to peak heights of about 500 relative fluorescent units (RFU)
- Purify the PCR products before starting the analysis

5.2 Setting up the Data Collection Software

Edit the run module as follows for the first run:

- In the **Module Manager** of the Data Collection Software click on **New** to open the **Run Module Editor** dialog box

Run Module 3kV_10s_450bp

| Parameter | Set up |
|-------------------------|----------------|
| Oven Temperature [°C] | <i>Default</i> |
| Poly Fill Volume | <i>Default</i> |
| Current Stability [µA] | <i>Default</i> |
| PreRun Voltage [kV] | <i>Default</i> |
| PreRun Time [s] | <i>Default</i> |
| Injection Voltage [kV] | 3.0 |
| Injection Time [s]* | 10 |
| Voltage Number of Steps | <i>Default</i> |
| Voltage Step Interval | <i>Default</i> |
| Data Delay Time [s] | <i>Default</i> |
| Run Voltage [kV] | <i>Default</i> |
| Run Time [s]** | 1500 |

* Deviating from the standard settings, the injection time may range between 1 and 20 s depending on the type of sample. If references samples with very high signal intensities are recorded, a shorter injection time may be selected in order to avoid pull-up peaks. For samples with low DNA content an injection time of up to 20 s may be necessary.

** Depending on the analysis conditions the run time for Mentype® **DIPscreen** was modified in order to be able to analyse fragments with lengths of up to **450 bp**.

- Click on **Save As**, enter the name of the new module (e.g. 3kV_10s_450bp) and confirm with **OK**
- Click on **Close** to exit the **Run Module Editor**

Starting the run

- Place the prepared 96-well plate on the autosampler tray
- In the **Protocol Manager** of the Data Collection software, click on **New** in the **Instrument Protocol** window to open the **Protocol Editor** dialog box

Instrument Protocol

| Protocol Editor | Set up |
|-----------------|----------------|
| Name | enter a name |
| Type | REGULAR |
| Run Module* | 3kV_10s_450bp |
| Dye Set | Any5Dye |

* parameter see above

- Click on **OK** to exit the **Protocol Editor**

Prior to each run, it is necessary to create a plate definition as follows:

- In the **Plate Manager** of the Data Collection software click on **New** to open the **New Plate Dialog** box

Plate Editor (I)

| New Plate Dialog | Set up |
|----------------------------|-------------------------------|
| Name | e.g. Plate_BT5_Date |
| Application | Select GeneMapper Application |
| Plate Type | 96-Well |
| Owner Name / Operator Name | ... |

- Click on **OK**. A new table in the **Plate Editor** will open automatically

Plate Editor (II)

| Parameter | Set up |
|-----------------------|---|
| Sample Name | Enter name for the samples |
| Priority | e.g. 100 (Default) |
| Sample Type | Sample or allelic ladder |
| Size Standard | e.g. SST-BTO_60-450bp |
| Panel | e.g. DIPscreen_Panels_v0 |
| Analysis Method | e.g. Analysis_DIPscreen_3130_200rfu |
| Snp Set | - |
| User-defined 1-3 | - |
| Results Group 1 | (select results group) |
| Instrument Protocol 1 | Run36_POP4_BT5_25min (setting described before) |

- Click into the column header to select the entire column, select **Edit** → **Fill Down** to apply the information to all selected samples and click on **OK**
- In the **Run Scheduler**, click on **Find All**, select **Link** to link the reaction plate on the autosampler up with the newly created plate record (position A or B) and start the run
- During the run, view **Error Status** in the **Event Log** or examine the quality of the raw data for each capillary in the **Capillaries Viewer** or the **Cap/Array Viewer**
- View data as overview in **Run History** or **Cap/Array Viewer** of the Data Collection software. Run data are saved in the **Run Folder** of the previously chosen **Result Group**

5.3 Analysis parameter / analysis method

The recommended analysis parameters are:

| | |
|---------------------------|--|
| Peak Detection Algorithm | Advanced |
| Ranges | Analysis: Full Range Sizing: All Sizes |
| Smoothing and Baselineing | Smoothing: Light Baseline Window: 51 pts |
| Size Calling Method | Local Southern Method |
| Peak Detection | Peak Amplitude Thresholds B:* Y:* G:* R:* O:* Min. Peak Half Width: 2 pts Polynomial Degree: 3 Peak Window Size: 15 pts Slope Thresholds: 0.0 |

Recommend settings in the worksheet **Allele** are:

Amelogenin Cutoff** 0.1

Recommend settings in the worksheet **Peak Quality** are:

| | |
|----------------------|----------------------------|
| Heterozygote balance | Min peak height ratio: 0.1 |
| Allele number | Max expected alleles: 2 |

* The peak amplitude threshold (Cutoff value) corresponds to the minimum peak height that will be detected from the GeneMapper® ID/ID-X software. The thresholds are usually 50-200 RFU and should be determined individually by the laboratory. Recommendation: The minimal peak height should be three times higher than the background noise of the baseline.

** All DIPs will be examined by GeneMapper® ID/ID-X software like Amelogenin.

6. Electrophoresis using the ABI PRISM® 3500/3500xL Genetic Analyzer

For detailed instructions on instrument setup, spectral calibration, or application of the Applied Biosystems 3500 Series Data Collection Software version 1.0 and the GeneMapper® ID-X software version 1.2, refer to the *Applied Biosystems 3500/3500xL Genetic Analyzers User Guide*.

The system with 8 capillaries is named AB 3500 and the system with 24 capillaries is named AB 3500xL.

The virtual **filter set Any5Dye** shall be used for the combined application of five fluorescent labels **6-FAM, BTG, BTY, BTR, and BTO** (the matrix standard will be called **BT5** hereinafter).

Material

| | |
|-----------|---|
| Capillary | 36 cm Capillary Array for 3500/3500xL |
| Polymer | POP-4® Polymer for 3500/3500xL |
| Buffer | 10x Genetic Analyzer Buffer with EDTA for 3500/3500xL |

6.1 Spectral calibration / matrix generation

Prior to conducting DNA fragment size analysis, it is necessary to perform a spectral calibration with the fluorescent labels **6-FAM, BTG, BTY, BTR, and BTO** for each analyzer. The calibration procedure creates a matrix that is used to correct the overlap of fluorescence emission spectra of the dyes.

Spectral calibration comprises the following steps:

- Preparation of spectral calibration standards
- Loading the standards to the multi-well reaction plate (one sample per capillary)
- Preparation of instrument and creating a Dye Set BT5
- Performing a spectral calibration run and checking the matrix

Setting up the spectral calibration standards

Example for 8 capillaries/ABI 3500

| Component | Volume |
|--|----------|
| Hi-Di™ Formamide | 108.0 µl |
| Matrix standard BT5 | 9.0 µl |
| - Load 12 µl of the mix to a 96-well reaction plate, e.g. position A1-H1 - Denaturation for 3 min at 95 °C - Cool down to 4 °C and place samples in the autosample tray | |

Example for 24 capillaries/ABI 3500xL

| Component | Volume |
|--|----------|
| Hi-Di™ Formamide | 300.0 µl |
| Matrix standard BT5 | 25.0 µl |
| - Load 12 µl of the mix to a 96-well reaction plate, e.g. position A1-H1, A2-H2 and A3-H3* - Denaturation for 3 min at 95 °C - Cool down to 4 °C and place samples in the autosample tray | |

* When using a 384-well plate, load 10 µl of the mixtures to columns 1, 3, and 5 in rows A, C, E, G, I, K, M, and O.

Performing a spectral calibration run

- Place the multi-well plate on the autosampler tray
- Now prepare the instrument and specific spectral calibration run settings

Preparation of the instrument

Before starting the spectral calibration process ensure that the spatial calibration has been performed. This process is necessary if a new capillary array was installed before and is described in detail in the *Applied Biosystems 3500/3500xL Genetic Analyzers User Guide*.

Preparation of dye set BT5

Prior to the spectral calibration, a dye set for the Matrix Standard BT5 needs to be setup.

1. To create a new dye set, go to **Library** and select **Analyze**, followed by **Dye Sets** and click **Create**.
2. Enter a **Dye Set Name**, e.g. BT5.
3. Select **Matrix Standard** as a chemistry and **AnyDye Template** as a Dye Set Template.
4. Disable **Purple** in the field **Arrange Dyes**. Ensure that all other colors are enabled.
5. Under **Calibration Peak Order** the colors need to be arranged as follows: **5 – blue, 4 – green, 3 – yellow, 2 – red, and 1 – orange**.
6. Do not alter the **Parameter** settings.
7. Click **Save** to confirm the changes.

Create New Dye Set

Setup a Dye Set

* Dye Set Name: BT5

* Chemistry: Matrix Standard

* Dye Set Template: AnyDye Template

Arrange Dyes

| | | | | | | |
|------------------------|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|
| Dye Selection | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> |
| Reduced Selection | <input checked="" type="checkbox"/> |
| Calibration Peak Order | 5 | 4 | 3 | 2 | 0 | 1 |

Parameters

The parameters will be used for instruments configured with 36cm capillary array and polymer POP4

Matrix Condition Number Upper Limit: 20.0

Locate Start Point: * After Scan: 300 * Before Scan: 5000

* Limit Scans To: 20000

Sensitivity: 0.1

* Minimum Quality Score: 0.8

Notes

Matrix Std. BT5 multi cap.

Close Save

Fig. 5 Setup for dye set BT5

- In the **Protocol Manager** of the Data Collection software click on **New** in **Instrument Protocol** to open the **Protocol Editor** dialog box

Performing a spectral calibration run

Once the multiwell plate containing the spectral calibration mixture is placed in the autosampler tray the spectral calibration process can be started.

1. To access the Spectral Calibration screen, select **Maintenance** on the Dashboard of the 3500 Series Data Collection software.
2. The number of wells in the spectral calibration plate and their location in the instrument must be specified.
3. Select **Matrix Standard** as a chemistry standard and **BT5** for dye set.
4. (Optional) Enable **Allow Borrowing**.
5. Click **Start Run**.




Fig. 6 Electropherogram of spectral calibration with matrix standard BT5 on an ABI 3500

Matrix check

- The quality value (**Q value**) of each capillary must be greater than 0.8 and the condition number range (**C value**) must be between 1 and 20
- Check the matrix samples for a flat baseline. As shown in the figure above, there should be five peaks with peak heights of about 1000-5000 RFU (Y-axis) in each matrix sample (optimal range: 2000-4000 RFU)
- A successful calibration will be displayed in green in **Overall** and for each capillary
- If all capillaries have passed the test, **Accept Results**
- If calibration failed, **Reject Results** and refer to **spectral calibration troubleshooting** of Applied Biosystems 3500/3500xL Genetic Analyzer User Guides

6.2 Sample preparation

| Component | Volume |
|--|---------|
| Hi-Di™ Formamide | 12.0 µl |
| DNA Size Standard 550 (BTO) | 0.5 µl |
| prepare 12 µl of the mix (formamide + DNA size standard) for all samples | |
| add 1 µl PCR product (diluted if necessary) or allelic ladder | |
| - Denaturation for 3 min at 95 °C | |
| - Cool down to 4 °C and place the samples on the autosample tray | |

Since injections take place simultaneously on all capillaries, 8 or 24 samples must be pipetted on the plate of multi-capillary analysers. If fewer samples are analysed empty positions need to be filled with 12 µl Hi-Di™ Formamide.

To ensure a reliable allelic assignment on multi-capillary analysers, several ladders should be run.

Room temperature may influence the performance of PCR products on multi-capillary instruments, so that shoulder peaks or split peaks occur especially at low

temperatures. Pay attention to keeping ambient conditions as recommended by the instrument manufacturer. Optimal will be a stable room temperature > 22 °C.

Signal intensities

Options to increase the signal intensity:

- Reduce the volume of the DNA Size Standard 550 (BTO) to peak heights of about 500 relative fluorescent units (RFU)
- Purify the PCR products before starting the analysis

6.3 Setting up a run

For the first run using the Mentype® **DIPscreen** Kit you will need to setup a number of protocols within the 3500 Series Data Collection Software.

Create Instrument protocol

- Go to **Library** and select **Analyze / Instrument protocol** and click **Create**
- Change the parameters according to the table below

Instrument protocol for Mentype® **DIPscreen**®

| Parameter | Set up |
|------------------------|------------------------|
| Application Type | HID / Microsatellite |
| Capillary Length | <i>Default</i> |
| Polymer | <i>Default</i> |
| Dye Set | BT5 |
| Run Module | <i>Default</i> |
| Protocol Name | e.g. Mentype DIPscreen |
| Oven Temperature [°C] | <i>Default</i> |
| Run Voltage [kV] | <i>Default</i> |
| Injection Voltage [kV] | 3.0 |
| Run Time [s]** | 1500 |
| PreRun Time [s] | <i>Default</i> |
| Injection Time [s]* | 10 |
| Data Delay Time [s] | <i>Default</i> |
| Advanced Options | <i>Default</i> |

* Deviating from the standard settings, the injection time may range between 1 and 20 s depending on the type of sample. If reference samples with very high signal intensities are recorded, a shorter injection time may be selected in order to avoid pull-up peaks. For samples with low DNA content an injection time of up to 20 s may be necessary.

** Depending on the analysis conditions, the run time for Mentype® **DIPscreen** was modified in order to analyse fragments with lengths of up to **450 bp**.

- Click on **Save** to confirm the settings

Create Size Standard

- Go to **Library** and select **Analyze / Size Standards** and click **Create**
- Change the parameters according to the table below

| Parameter | Set up |
|---------------|---------|
| Size Standard | BTO_550 |
| Dye Color | Orange |

The DNA Size Standard 550 (BTO) should be used with the following lengths of fragments: **60, 80, 90, 100, 120, 140, 160, 180, 200, 220, 240, 250, 260, 280, 300, 320, 340, 360, 380, 400, 425, 450, 475, 500, 525, and 550** bp.

- Click on **Save** to confirm the settings

Create QC (Size Calling) Protocol

- Go to **Library** and select **Analyze / QC (Size Calling)** and click **Create**
- Change the parameters according to the table below

| Parameter | Set up |
|---------------|----------------------|
| Protocol Name | enter a name |
| Size Standard | BTO_550 (from above) |
| Sizercaller | Size Caller v.1.1.0 |

- Go to **Analysis Settings / Peak Amplitude Threshold** and **disable purple**.
All other colours should be enabled.
- Keep all other settings as Default
- Click on **Save** to confirm the settings

Create an Assay

- Go to **Library** and select **Manage / Assays** and click **Create**
- Change the parameters according to the table below

| Parameter | Set up |
|---------------------|------------------------|
| Assay Name | e.g. Mentype DIPscreen |
| Color | Default |
| Application Type | HID |
| Instrument Protocol | e.g. Mentype DIPscreen |
| QC Protocols | e.g. BTO_550 |

- Click on **Save** to confirm the settings

Starting the run

- Place the prepared multi-well plate on the autosampler tray
- In the **Dashboard** of the Data Collection Software, click **Create New Plate**
- Go to **Define Plate Properties** and select **Plate Details**
- Change the parameters according to the table below

Plate Details

| Property | Set up |
|------------------|--------------|
| Name | enter a name |
| Number of Wells | 96 or 384 |
| Plate Type* | HID |
| Capillary Length | 36cm |
| Polymer | POP4 |

- Click **Assign Plate Contents** to confirm the settings
- Define well position of each sample or ladder for data collection and processing by entering sample names
- Assign an **Assay** (required) a **File Name Conventions** and a **Result Group** to all named wells in the plate
- Click **Link the plate for Run** and enter Run Name
- Click **Start Run**

7. Analysis

For general instructions on automatic sample analysis, refer to the *GeneScan®* or *GeneMapper® ID* or *GeneMapper® ID-X Software User's Manual*.

Note: Within the Mentype® DIPscreen the red panel should be faded out.

Finding the exact lengths of the amplified products depends on the device type, the conditions of electrophoresis, as well as the DNA size standard used. The DNA Size Standard 550 (BTO) shall thus be used with the following lengths of fragments: **60, 80, 90, 100, 120, 140, 160, 180, 200, 220, 240, 250, 260, 280, 300, 320, 340, 360, 380, 400, 425, 450, 475, 500, 525, and 550** bp.




Fig. 7 Electropherogram of the DNA Size Standard 550 (BTO), fragments with lengths in bp

Note: The provided template files for the DNA size standard SST-BTO_60-450bp can be applied for the evaluation and analysis of the Mentype® **DIPscreen** using the *GeneMapper® ID* or *ID-X Software*.

7.1 Biotype® template files

Allele allocation should be carried out with suitable analysis software, e.g. the GeneMapper® ID/ID-X software in combination with the Mentype® DIPscreen template files from Biotype or the Chimeris™ **Monitor** Software of Biotype Diagnostic GmbH. Template files are available from our homepage or as CD-ROM on request.

Recommended Biotype® templates for GeneMapper® ID/ID-X Software are:

| | | |
|-----------------|--|--------------------|
| Panels | DIPscreen_Panels_v0/v0X | or higher versions |
| BinSets | DIPscreen_Bins_v0/v0X | or higher versions |
| Size Standard | SST-BT0_60-450bp | |
| Analysis Method | Analysis_DIPscreen_310_200rfu Analysis_DIPscreen_310_1000rfu Analysis_DIPscreen_3130_200rfu Analysis_DIPscreen_3130_1000rfu | |
| Plot Settings | PlotsBT5_4dyes | |
| Table Settings | Table for 2 alleles | |

Panels and BinSets always have to be used whereas the other template files are optional.

The prepared Biotype® templates for GeneMapper® ID/ID-X Software were generated for POP4® runs. In case of using other polymer types changes may necessary on Panels and Bins or within the Analysis Method before analyzing the data.

For detailed instruction please refer to the instruction Biotype® Template Files für GeneMapper® which can be downloaded from our homepage (www.biotype.de).

Important Note: Import and allele calling with provided template files is only guaranteed using GeneMapper® ID/ID-X software. If GeneMapper® software is applied you may experience import problems using some template files. You may have to adjust Panels and Bins with one ore more runs of the allelic ladder on your specific instrument setup. Contact us for support (support@biotype.de).

General procedure for the analysis

1. Check the DNA size standard
2. Check the allelic ladder
3. Check the positive control
4. Check the negative control
5. Analyse and interpret the sample data

7.2 Controls

The Control DNA XY13 of the test kit and other commercially available DNA from standard cell lines represent the following alleles:

Table 2. Allele determinations of Mentype® DIPscreen

| Locus | Control-DNA XY13 | ATCC K-562 | CCR 9947A | CCR 9948 | CCR 3657 |
|--------|------------------|------------|-----------|----------|----------|
| AM | XY | XX | XX | XY | XY |
| HLD106 | +/+ | -/- | +/+ | +/+ | +/+ |
| HLD70 | -/+ | -/+ | +/+ | -/+ | -/- |
| HLD84 | -/+ | +/+ | -/- | -/+ | -/- |
| HLD103 | +/+ | -/- | -/+ | +/+ | -/+ |
| HLD104 | -/+ | -/- | -/+ | +/+ | -/- |
| HLD116 | -/+ | +/+ | -/- | -/+ | -/- |
| HLD112 | -/+ | +/+ | -/+ | -/+ | -/+ |
| HLD307 | +/+ | +/+ | -/+ | +/+ | +/+ |
| HLD310 | +/+ | -/+ | -/+ | -/- | -/+ |
| HLD110 | -/+ | -/+ | -/+ | -/+ | -/+ |
| HLD133 | -/+ | -/- | +/+ | +/+ | -/+ |
| HLD79 | +/+ | +/+ | +/+ | -/+ | +/+ |
| HLD105 | -/+ | -/- | -/+ | -/+ | -/+ |
| HLD140 | +/+ | +/+ | -/- | -/+ | +/+ |
| HLD163 | +/+ | -/+ | -/+ | +/+ | -/+ |
| HLD91 | -/+ | -/+ | -/- | -/- | -/+ |
| HLD23 | -/+ | +/+ | -/- | -/+ | -/+ |
| HLD88 | +/+ | -/- | -/- | -/+ | +/+ |
| HLD101 | -/+ | -/+ | -/+ | -/+ | -/+ |
| HLD67 | -/+ | -/+ | +/+ | +/+ | +/+ |
| HLD301 | -/+ | -/+ | -/+ | -/+ | -/- |
| HLD53 | +/+ | -/- | -/+ | +/+ | -/- |
| HLD97 | -/- | -/- | -/+ | -/+ | +/+ |
| HLD152 | -/- | +/+ | +/+ | -/+ | +/+ |
| HLD128 | -/+ | -/+ | -/+ | -/- | -/+ |
| HLD134 | -/+ | -/- | +/+ | +/+ | -/- |
| HLD305 | -/+ | -/+ | -/+ | +/+ | -/+ |
| HLD48 | -/+ | +/+ | +/+ | -/+ | +/+ |
| HLD114 | +/+ | -/- | -/- | +/+ | -/+ |
| HLD304 | +/+ | -/- | -/+ | -/+ | -/- |
| HLD131 | +/+ | -/+ | -/- | -/+ | +/+ |
| HLD38 | +/+ | -/+ | -/+ | +/+ | +/+ |
| HLD82 | +/+ | +/+ | +/+ | -/+ | +/+ |

The reference DNA K-562 is available from ATCC. DNA 9947A, 9948 and 3657 are available from Coriell Cell Repositories.

7.3 Lengths of fragments and alleles

Table 3 show the fragment lengths of individual alleles that refer to the DNA Size Standard 550 (BTO). All analyses have been performed on an ABI PRISM® 3130 Genetic Analyzer with POP4® polymer. Different analysis instruments, DNA size standards or polymers may result in different fragment lengths. In addition, a visual alignment with the allelic ladder is recommended.

Scaling

Horizontal: 70-430bp (see Fig. 8 and 9)

Vertical: Depending on signal intensity

Figure 9




Fig. 9 Electropherogram of the allelic ladder Mentype® DIPscreen. Analysis was performed on an ABI PRISM® 3130 Genetic Analyzer with the DNA Size Standard 550 (BT). Allele assignment was performed using the GeneMapper® ID Software and the Mentype® DIPscreen template file.

Table 3. Fragment lengths of the Mentype® DIPscreen allelic ladder analysed on an ABI PRISM® 3130 Genetic Analyzer with POP4® (FAM, BTG, BTY panel)

| Marker/FAM | -DIP [bp]* | +DIP [bp]* | Marker/BTG | -DIP [bp]* | +DIP [bp]* |
|------------|------------|------------|------------|------------|------------|
| AM | 77 (X) | 80 (Y) | HLD91 | 84 | 88 |
| HLD106 | 91 | 98 | HLD23 | 107 | 113 |
| HLD70 | 104 | 108 | HLD88 | 118 | 128 |
| HLD84 | 112 | 117 | HLD101 | 131 | 135 |
| HLD103 | 129 | 138 | HLD67 | 140 | 148 |
| HLD104 | 153 | 1157 | HLD301 | 172 | 182 |
| HLD116 | 170 | 192 | HLD53 | 190 | 194 |
| HLD112 | 199 | 204 | HLD97 | 214 | 228 |
| HLD307 | 228 | 236 | HLD152 | 239 | 250 |
| HLD310 | 248 | 257 | HLD128 | 258 | 266 |
| HLD110 | 264 | 270 | HLD134 | 296 | 312 |
| HLD133 | 278 | 288 | HLD305 | 375 | 401 |
| HLD79 | 294 | 299 | | | |
| HLD105 | 302 | 310 | Marker/BTY | -DIP [bp]* | +DIP [bp]* |
| HLD140 | 318 | 333 | HLD48 | 78 | 83 |
| HLD163 | 344 | 358 | HLD114 | 159 | 177 |
| | | | HLD304 | 184 | 203 |
| | | | HLD131 | 208 | 220 |
| | | | HLD38 | 234 | 240 |
| | | | HLD82 | 314 | 352 |

* rounded to integer

8. Interpretation of results

As mentioned above, post PCR analysis and automatic allele allocation with suitable analysis software ensure a precise and reliable discrimination of alleles.

Pull-up peaks

Pull-up peaks may occur if peak heights are outside the linear detection range (>3000 RFU), or if an incorrect matrix was applied. They appear at positions of specific peaks in other colour channels, typically with lower signal intensities. Peak heights should not exceed 3000 RFU in order to prevent pull-up peaks.

Template-independent addition of nucleotides

Because of its terminal transferase activity, the Multi Taq DNA Polymerase tends to add an adenosine radical at the 3'-end of the amplified DNA fragments. The artefact peak is one base shorter than expected (-1 bp peaks). All Biotype® primers are designed to minimise these artefacts. Artefact formation is further reduced by the final extension step of the PCR protocol at 68 °C for 60 min. Peak height of the artefact correlates with the amount of DNA. Laboratories should define their individual limits for analysis of the peaks.

Artefacts

Room temperature may influence the performance of PCR products on multi-capillary instruments, shoulder peaks or split peaks occur. Furthermore, automated assignment could be influenced in some cases. If these effects occur we recommend injecting the sample again at higher room temperature and maybe using more than one allelic ladder sample per run.

Influence of polymers

Mentype® **DIPscreen** was validated and certified for the analysis on POP4® polymer. The use of other polymers (e.g. POP7™ or POP6™) might influence the run behaviour of specific PCR products. Furthermore background noise might increase through different behaviour of free fluorescent dyes.

9. References

- Alizadeh M, Bernard M, Danic B, Dauriac C, Birebent B, Lapart C, Lamy T, Le Prise PY, Beauplet A, Bories D, Semana G, Quelvennec E. (2002)** Quantitative assessment of hematopoietic chimerism after bone marrow transplantation by real-time quantitative polymerase chain reaction. *Blood* 99, 4618-4625.
- Chen DP, Tseng CP, Wang WT, Wang MC, Tsai SH, Sun CF (2011)** Real-time biallelic polymorphism-polymerase chain reaction for chimerism monitoring of hematopoietic stem cell transplantation relapsed patients. *Clin Chim. Acta* 412, 625-630.
- Harries LW, Wickham CL, Evans JC, Rule SA, Joyner MV, Ellard S (2005)** Analysis of haematopoietic chimerism by quantitative real-time polymerase chain reaction. *Bone Marrow Transplant.* 35, 283-290.
- Masmas TN, Madsen HO, Petersen SL, Ryder LP, Svejgaard A, Alizadeh M, Vindelov LL (2005)** Evaluation and automation of hematopoietic chimerism analysis based on real-time quantitative polymerase chain reaction. *Biol Blood Marrow Transplant.* 11, 558-566.
- Mills RE, Luttig CT, Larkins CE, Beauchamp A, Tsui C, Pittard WS, Devine SE (2006)** An initial map of insertion and deletion (INDEL) variation in the human genome. *Genome Res* 16 (9):1182-1190, 2006.
- Qin XY, Li GX, Qin YZ, Wang Y, Wang FR, Liu DH, Xu LP, Chen H, Han W, Wang JZ, Zhang XH, Li JL, Li LD, Liu KY, Huang XJ (2011)** Quantitative assessment of hematopoietic chimerism by quantitative real-time polymerase chain reaction of sequence polymorphism systems after hematopoietic stem cell transplantation. *Chin Med J (Engl.)* 124, 2301-2308.
- Weber JL, David D, Heil J, Fan Y, Zhao C, Marth G (2002)** Human diallelic insertion/deletion polymorphisms. *Am J Hum Genet* 71(4):854-862.
- Wilhelm J, Reuter H, Tews B, Pingoud A, Hahn M (2002)** Detection and quantification of insertion/deletion variations by allele-specific real-time PCR: application for genotyping and chimerism analysis. *Biol Chem* 383, 1423-1433.

10. Explanation of Symbols




Manufacturer



Date of manufacture



Batch code



**Contains sufficient reagents for
<N> tests**



**Consult instructions (handbook)
for use**



Use by



Temperature limitations



Catalogue number



In-Vitro-Diagnostics

Specification Mentype® DIPscreen PCR-Amplification kit

A Analytical Validation

A a) Determination of the Standard Reaction and batch-specific Tolerance

Objective: The standard reaction and the batch-specific tolerances with respect to the absolute signal heights (RFU), the balance of the signal heights of the multiplex PCR and the baseline were determined.

Methodology: The test kit contains the control DNA XY13 of a healthy donor, which is heterozygous in 17 DIP systems and amelogenin (AM). The standard reaction (28 PCR cycles) was carried out with this control DNA in the nominal concentration of 1 ng in quadruple determination. Four blank values (no template control, NTC) without DNA were also carried out.

Results: For the batch-specific mixing of the PCR primers, the following specifications were established: Signal heights of 1000-5000 RFU were obtained using an ABI PRISM® 3130 Genetic Analyzer. The fluctuations for signal heights of heterozygous systems allowed a maximum of 50 % of the guide value. In the scaling range, no unspecific signals 200 RFU were determined (baseline) for the blank values.

A b) Genotyping Accuracy

Objective: The accuracy of the allele assignment should be statistically secured under standard conditions. The testing examined the automatic allele calling with the allelic leader. Further the concordance of the allele assignment compared to the pre-typing of test DNAs by other methods (other PCR kits, Direct sequencing, etc.) using the GeneMapper ID software. Based on the results, the test-specific device settings for the genotyping by means of capillary electrophoresis (bins and panels) for the analysis templates of the DNA sequencer are defined.

Methodology: 100 pre-typed human DNAs from donors (whole blood, cheek swabs) were investigated in single determination. In addition, a blank was carried without DNA. Acceptance criterion was defined as full profiles with peak heights ≥ 200 RFU (manual evaluation) [3; 4].

Results: After determination of the test-specific device settings the correct genotype was assigned to all DNA samples for all HLD systems and the amelogenin marker.

A c) Analytical Specificity

Objective: The investigations served to exclude false-positive results due to cross-reactivity with selected non-human DNA samples. In clinical practice, however, non-human DNA can be largely excluded due to sterile sampling.

Method: 2.5 ng genomic DNA from *Bos Taurus* (Cattle), *Sus scrofa domestica* (Pig), *Canis lupus familiaris* (dog), *Felis catus* (cat) and *Oryctolagus cuniculus* (domestic rabbit) were tested. The DNA from animals was derived from blood samples, which were provided as a residual material of veterinary studies.

Results: No cross-reactivity was detected in the allele area (<200 RFU).

A d) Analytical Sensitivity

Objective: The investigations were used to determine the analytical detection limit (sensitivity).

Method: A dilution series with 1 ng to 65 pg of reference DNA was tested in quadruplicate. As an acceptance criterion, complete DNA profiles with ≥ 100 RFU were defined.

Results: A detection limit of 200 pg of genomic DNA was determined.

A e) Assays Performance with Different PCR-Thermocyclers

Objective: PCR thermal cyclers of different manufacturers differ in their specifications. In particular, different heating and cooling rates as well as different temperature control techniques (hysteretic versus symmetrical settling of the actual values by the setpoint) can be observed.

Method: Testing of the standard reaction with control DNA in the nominal concentration of 1 ng was carried out with all thermocyclers (described below) in quadruplicate determinations with the same master mix. In addition, 2 blank samples without DNA were examined.

GeneAmp 9700 with Silverblock (Applied Biosystems®, Life Technology GmbH, Darmstadt), GeneAmp 9700 with Alu block (Life Technology GmbH, Darmstadt) and Eppendorf Mastercycler ep-S (Eppendorf AG, Hamburg)

Results: No unspecific by-products ≥ 200 RFU were detected in the allele region. The deviation of the mean peak heights compared to the standard reaction was a maximum of 20% at a defined ramp of ≥ 2 ° C./sec.

A f) Mixed DNA Samples

Objective: The aim of chimerism analysis after allogeneic blood stem cell transplantation is the separate detection and the relative quantification of donor and recipient DNA. In order to detect the minimal residual disease, the smallest possible amounts of recipient DNA should be detected in the mixture. Therefore, different mixtures of two defined DNAs with different genotypes were produced in analytical validation.

Method: 10 independent mixtures of each two unrelated DNAs were prepared using the deficient DNA as 0%, 1%, 5%, 10%, 30%, 50% and 70%. Between two DNAs in the mixtures, an average of 13 DIP-Loci (12.8 ± 2.22) with informative alleles were added for evaluation. In each case 2 ng of the DNA mixtures was tested in the standard reaction. Signal heights of at least 50 RFU were evaluated.

Results: A detection limit of 1% could be achieved for the deficient DNA. This corresponds to the values 1 - 5%, which were achieved with forensic STR kits in the chimerism analysis [8-11].

A g) PCR Annealing Temperatures

Objective: To determine the robustness of the PCRs, temperature fluctuations for the primer attachment step (annealing) of the multiplex PCR were simulated. This temperature step is critical for the sensitivity and specificity of the PCR.

Method: The kit-specific annealing temperature of 60°C was varied by $\pm 1^\circ\text{C}$ and $\pm 2^\circ\text{C}$ in the standard reaction setting with control DNA and a nominal concentration of 1 ng. A triple determination with the same master mix performed.

Results: No nonspecific by-products ≥ 200 RFU were detected for $\pm 1^\circ\text{C}$. The average peak heights deviated from the standard reaction at $\pm 1^\circ\text{C}$ to a maximum $\pm 30\%$. At $+2^\circ\text{C}$ some systems (HLD 84, 103, 116, 112, 133, 105, 40, 67, 48), showed a reduction in performance. On system showed a complete profile loss (HLD 91) at that temperature.

A h) Fluctuation of PCR buffer Batches

Objective: The concentration ratios of the contents of the PCR buffer Reaction mix A (dNTPs, ion concentrations, in particular Mg^{2+}) are critical for sensitivity, specificity and the balance of the signals in multiplex PCRs. Therefore the robustness of the test is tested against batch fluctuations of the supplied PCR buffer.

Methods: 4 independent Reaction mix A batches were tested for performance in the standard reaction with control DNA of the nominal concentration of 1 ng.

Results: No unspecific by-products ≥ 200 RFU were detected. The deviation of averaged peak heights compared to the standard reaction was a maximum of 20%.

A i) In-use Stability

Objective: The stability of the reagents of the PCR kit was tested after repeated freezing and thawing.

Methodology: The kit reagents were subjected to a 20-fold freezing and thawing cycle. The freezing was carried out for at least 1 h at -20 °C. The mixture was thawed at room temperature and the reagents were homogenized by shaking before use. Subsequently, a standard reaction with control DNA of the nominal concentration of 1 ng and additional blank values without DNA was performed in triplicate determinations. The evaluation was carried out in comparison to a standard reaction without freezing and thawing cycle.

Results: The deviation of averaged peak heights compared to the standard reaction was a maximum of 20 % (in particular signal loss). Additional peaks >200 RFU were observed in the blank samples; however, no peak occurred in the allelic range of the kit (free fluorescence dyes in the BTG Panel).

B Clinical Performance Data

B a) Study Design, ethics and regulatory aspects

A clinical performance study was conducted according to §§ 20 - 24 of German Medizinproduktegesetz. The protocol was approved by the National Competent Authority BfArM according to § 7 German Verordnung über klinische Prüfungen von Medizinprodukten and by the institutional ethics committee. All participants gave written informed consent.

B b) Reference Methods

The performance of the kit was compared to the Short-Tandem-Repeats (STRs) based CE-IVD Mentype® Chimera PCR-Amplification kit (Biotype Diagnostic GmbH, Dresden, DE) [12]. Further, a cytogenetic differentiation of donor and recipient leukocytes by means of fluorescence in situ hybridization (FISH) was performed using the sex chromosome-specific CE-IVD CEP® X SpectrumOrange™ / Y SpectrumGreen™ Direct Labeled Fluorescent DNA Probe Kit (Abbott GmbH & Co KG, Wiesbaden, DE; usage in according to the manufacturer's data) [11].

B c) DNA-Extraction and Purification



DNA extraction from heparinized whole blood samples was performed with the QIAamp® DNA Blood Mini Kit (Qiagen GmbH, Hilden, DE), according to the manufacturer.

B d) Results

A total of 98 data sets of adult patients were collected on various days after allogeneic blood stem cell or bone marrow transplantation. The donor-recipient pairs differed in the genetic gender and were thus suitable for the sex chromosome-specific FISH [10]. Per PCR, at least 1.5 ng genomic DNA was used. Firstly, all informative STR or DIP systems of the donor-recipient pairs were determined and sex was confirmed by genotyping the amelogenin marker, which is part of the multiplex PCR. For the PCR results mean values of the signal heights of all informative STR- or DIP systems were used. The results of the concordance analysis are summarized in Fig.1.

Compared to cytogenetics, 11 samples analyzed with Mentype® **DIPscreen** showed a deviation of the donor portion of more than 5% (absolute) (see Fig. 1A). In 5 of these samples, cell numbers significantly less than 200 were counted for cytogenetics. However, according to the recommendations of the manufacturer of the FISH kit, at least 200 cells should be counted. According to practical recommendations, higher absolute cell numbers (500-1000) yield better cytogenetic results [11, 13]. In contrast to cytogenetics, the differences from Mentype® **DIPscreen** to the STRs-based multiplex PCR kit Mentype® **Chimera** were 7.9% (see Fig. 1B). Only 3 of the 98 measurement data sets showed a deviation of more than 5%.

Fig. 1 Concordance analysis of multiplex PCR Mentype DIPscreen in comparison to cytology (A) and multiplex PCR Chimera (B)



B e) References

- 1) **Wenz H, Robertson JM, Menchen S, Oaks F, Demorest DM, Scheibler D, Rosenblum BB, Wike C, Gilbert DA, Efcavitch JW.** High-precision genotyping by denaturing capillary electrophoresis. *Genome Res* 1998; 8: 69–80.
- 2) **Sgueglia JB, Geiger S, Davis J.** Precision studies using the ABI prism 3100 genetic analyzer for forensic DNA analysis. *Anal Bioanal Chem* 2003; 376: 1247–54.
- 3) **Gilder JR, Doom TE, Inman K, Krane DE.** Run-specific limits of detection and quantitation for STR-based DNA testing. *J Forensic Sci* 2007; 52: 97–101.
- 4) **Schneider PM, Fimmers R, Keil W, Molsberger G, Patzelt D, Pflug W, Rothämel T, Schmitter H, Schneider H, Brinkmann B.** The German Stain Commission: recommendations for the interpretation of mixed stains. *Int J Legal Med.* 2009; 123: 1–5.
- 5) **Wheeler DL, Barrett T, Benson DA, Bryant SH, Canese K, Chetvermin V, Church DM, DiCuccio M, Edgar R, Federhen S, Geer LY, Kapustin Y, Khovayko O, Landsman D, Lipman DJ, Madden TL, Maglott DR, Ostell J, Miller V, Pruitt KD, Schuler GD, Sequeira E, Sherry ST, Sirotkin K, Souvorov A, Starchenko G, Tatusov RL, Tatusova TA, Wagner L, Yaschenko E.** Database resources of the National Center for Biotechnology Information. *Nucleic Acids Res* 2007; 35 (Database issue): D5–12.
- 6) **Haas-Rochholz H, Weiler G.** Additional primer sets for an amelogenin gene PCR-based DNA–sex test. *Int J Legal Med* 1997; 110: 312–5.
- 7) **Peakall R and Smouse PE.** GenAlEx 6.5: Genetic analysis in Excel. Population genetic software for teaching and research—an update. *Bioinformatics* 2012; 28: 2537–2539.
- 8) **Thiede C, Lion T.** Quantitative analysis of chimerism after allogeneic stem cell transplantation using multiplex PCR amplification of short tandem repeat markers and fluorescence detection. Appendix: Method in focus. *Leukemia* 2001; 15: 303–6.
- 9) **Thiede C.** Diagnostic chimerism analysis after allogeneic stem cell transplantation: new methods and markers. *Am J Pharmacogenomics* 2004; 4: 177–87.
- 10) **Thiede C, Lion T.** Quantitative analysis of chimerism after allogeneic stem cell transplantation using multiplex PCR amplification of short tandem repeat markers and fluorescence detection. Appendix: Method in focus. *Leukemia* 2001; 15: 303–6.
- 11) **Buño I, Nava P, Simón A, González-Rivera M, Jiménez JL, Balsalobre P, Serrano D, Carrión R, Gómez-Pineda A, Díez-Martín JL.** A comparison of fluorescent in situ hybridization and multiplex short tandem repeat polymerase

chain reaction for quantifying chimerism after stem cell transplantation. *Haematologica* 2005; 90: 1373–9.

- 12) **Henke L, Muche M, Blaauw A, Van Eede PH, Martin W, Helmken C, Budowle B, Henke J.** Validation of a "new" short tandem repeat (STR) fluorescent multiplex system and report of population genetic data. *Clin Lab* 2007; 53:477–82.
- 13) **Mohr B, Koch R, Thiede C, Kroschinsky F, Ehninger G, Bornhäuser M.** CD34+ cell dose, conditioning regimen and prior chemotherapy: factors with significant impact on the early kinetics of donor chimerism after allogeneic hematopoietic cell transplantation. *Bone Marrow Transplant* 2004; 34: 949–54.

Notes

Notes