

# PROTRANS Cyclerplate System

## Protrans Cyclercheck

**REF 200 100, 200 105**

### Instruction

- Take the Protrans PCR workstation out of the freezer (-20°C).  
Take the Protrans **Buffer D (blue cap)** and **Protrans Buffer E (purple cap)** and **Protrans Buffer Y** out of the freezer (-20°C). Place the Protrans Buffer D, E in the PCR workstation.  
Thaw Protrans Buffer Y and place it in the PCR workstation.

Protrans SSP Cycler Program			
Initial denaturation	94°C	2 min	Hold
Denaturation	94°C	15 sec	10 cycles
Annealing and Extension	65°C	60 sec	
Denaturation	94°C	15 sec	20 cycles
Annealing	61°C	50 sec	
Extension	72°C	30 sec	
Hold	4°C	15 min	Hold
Ramp rate 1°C/sec.			

- Take the **Protrans Cyclercheck Cyclerplate** out of the refrigerator (2-8°C) and place it in the PCR workstation. If you want, you can cut the Protrans 96-well Cyclercheck Cyclerplate with the Protrans cutter in 2x 48-well or 4x 24-well Cyclerplates.
- Preparation of the **Master Mix**.  
For each 96-well Protrans Cyclercheck Cyclerplate pipette in a 1.5ml reaction tube:

Protrans Cyclercheck	Buffer D	Buffer Y	Taq Polymerase 5U/µl	Buffer E
	<b>276 µl</b>	<b>710 µl</b>	<b>6.5 µl</b>	<b>50 µl</b>

- Vortex the **Master Mix** very thoroughly and spin the tube shortly.
- Remove the label from the Protrans Cyclercheck Cyclerplate and mark the plate.
- Dispense of the Master Mix** with a stepper:  
**10µl in all 96 positions** of the Protrans **Cyclercheck Cyclerplate**.
- Shake down the red Master Mix drops in the Protrans Cyclercheck Cyclerplate.
- Close the Protrans Cyclercheck Cyclerplate with the Protrans 96-well coverplate or strips and seal it very carefully by pressing the caps into the wells.
- Place Protrans Cyclercheck Cyclerplate directly in the thermocycler, final volume 10µl, and start the amplification immediately or store the covered cyclerplate at 2-8°C and start the amplification within 2 hours.

